



1)




Select the “Gear” icon on the side navigation bar to display order set screen. Select the desired lab for the order set, then select “Create Orderset.”

2)




Enter the desired name for the order set and then select “Start Adding Tests.”

3)




Select the desired test and move them into the orderset via the “Move” button highlighted by #2 arrow. Once complete select “Save Orderset.”

4)



The newly created order set will now show in your settings screen. Ordersets that are labeled “Private’ can only be utilized by you.

5)



The created orderset will now display on the lab order screen near the top of the corresponding labs test compendium.

LabSavvy IQ Support & Training Resources:



Operationsupport@labsavvy.com



[LabSavvy Training & Document Resources](#)



[Schedule a Training Session](#)



[Video Tutorial](#)