



1)




Select the “Patient Icon” on the side navigation panel to open the patient database.

2)




Select “New Patient” in the top right of the screen.

3)




Complete the indicated required information for the patient to be registered.

4)



Select “Save Changes” on the bottom right of the screen to save the patient profile & information.

5)



The patient is now available in your patient database and lab orders can be placed.

LabSavvy IQ Support & Training Resources:



Operationsupport@labsavvy.com



[LabSavvy Training & Document Resources](#)



[Schedule a Training Session](#)



[Video Tutorial](#)